


# CircAid<sup>®</sup> Systems

## T-3 M<sup>™</sup> Sizing Chart



T-3 M <sup>™</sup> Sizing Chart			
Measure ankle at narrowest point, calf at widest. If both measurements do NOT fall within one of these size ranges, choose size by the patient's calf measurement. Always use the undersock or stockinette when testing a size.			
LENGTH	SIZE	ANKLE	CALF
11" & 13"	Small	7½ - 8½" (19-21.5cm)	10 - 14" (25-35.5 cm)
11" & 13"	Medium	8½ - 10½" (21.5-26.5cm)	13 - 17" (33-43 cm)
11" & 13"	Large	10-12½" (25-32 cm)	16 - 20" (40.5-51 cm)
13"	X-Large	12-14½" (30.5 - 37 cm)	19 - 23" (48 - 58.5cm)

1 800 CIRCAID **CircAid<sup>®</sup>** 1-800-247-2243  
**MEDICAL + PRODUCTS**